

CR / DR Cassette Cleaning Log

_						
1)	Δ'	١,	11	~	$^{\circ}$	•
D	C	v	ı۷	٠,	C	•

For CR panels, inspect and clean both the outside cassette casings and the interior image phosphors quarterly. Clean the cassette image phosphors with manufacturer recommended cleaner. For DR panels, wipe down the panel casing only.

Date (dd/mm/yr)	Test	Pass (x)	Fail (x)	Comments	Initials
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				
	Cassettes Cleaned				